

**Product Expiration Date Policy**

Current lots are identified on the ChromaDex™ Extranet accessed through [www.chromadex.com](http://www.chromadex.com). In many cases, the previous lot may still be considered valid after the new lot is introduced. Ordinarily, a previous lot is considered as valid for about one year after the current lot enters distribution. It is the responsibility of each user to determine that their lot is suitable when used.

ChromaDex™ will not ship items with an expiration date of less than 6 months leading up to the listed expiration date on the certificate of analysis, unless pre-approved by the customer.

When standards approach the lot's expiration date and inventory still exists, they will be re-evaluated by the ChromaDex™ Analytics Laboratory.

- The standard will be re-analyzed by HPLC/GC (or whichever analysis is applicable).
- A comparison will be made to the data listed on the certificate of analysis.
- When the re-assay results fall within 5% of the listed HPLC/GC purity, a reverified date will be noted on the CofA and the HPLC/GC purity data on the certificate of analysis will not be revised.
- When the re-assay results fall outside of the 5% HPLC/GC purity, the lot is removed from inventory and the expiration date is not extended.
- RG lots are re-evaluated by MS or NMR verification.

Standards may be re-evaluated prior to the 6 month window leading to expiration.

**Contact your ChromaDex™ Sales Representative  
for more information on pricing and timing of  
custom lot reverification.**

